
Patient Information

AESCULAP®
Silver Graft

B BRAUN
SHARING EXPERTISE

AESCULAP® Silver Graft
Patient information

en	Patient information	3
de	Patienteninformation	4
fr	Information du patient	5
es	Información para el paciente	6
it	Informazioni per il paziente	7
pt	Informações do doente	8
nl	Informatie ten behoeve van de patiënt	9
da	Patientinformation	10
no	Pasientinformasjon	11
sv	Patientinformation	12
fi	Tiedon välittäminen potilaalle	13
et	Patsiendi teavitamine	14
lv	Informācija pacientam	15
lt	Pacientų švietimas	16
cs	Informace pro pacienty	17
pl	Informacja dla pacjenta	18
sk	Informácia pre pacientov	19
hu	A betegek oktatása	20
sl	Izobraževanje bolnikov	21
hr	Edukacija pacijenata	22
ro	Instruirea pacientului	23
bg	Образование на пациента	24
tr	Hasta bilgisi	25
el	Εκπαίδευση ασθενή	26

Silver Graft

AESCULAP® Patient Information

en

What does your body come into contact with after Silver Graft implantation?

Materials and substances


Silver Graft is a warp-knitted, double-velour vascular prosthesis, made of polyethylene terephthalate (polyester), for reconstructive procedures on the aorta and arteries in the peripheral area. Silver Graft is impregnated with absorbable modified bovine gelatine for temporary sealing. The outer surface of Silver Graft is coated with a metallic silver layer to inhibit microbial adhesion. Other ingredients are glycerol and hexamethylene diisocyanate (crosslinker for gelatine), Palladium and Titanium. In case of Silver Graft Helix additionally a polypropylene thread is included to protect the prosthesis against mechanical pressure from outside and increase the kink-stability.

How long will your implanted Silver Graft last?

Expected lifespan

The use of gelatine-coated, silver-impregnated polyester grafts such as Silver Graft as effective arterial substitutes has been assessed in both in vivo (animal) and in vitro (laboratory) experiments (1-4). Polyester is commonly used for vascular grafts in procedures such as bypass surgeries (5). Its durability and compatibility with the human body make it a preferred choice for many surgeons. The implant is not broken down by the body.

Therefore, Silver Graft remains implanted for the rest of the patient's life, if no adverse events occur. After implantation, the body's own connective tissue cells migrate into the implant.

What you and your doctor should look out for!

Interactions with environmental conditions

Magnetic resonance imaging (MRI) can be safely done in patients with implanted Silver Graft vascular prostheses under the following conditions:

- Static magnetic field: 1.5 T and 3 T
- Maximum spatial field gradient: 6.6 T/m
- Whole body averaged specific absorption rate (SAR): 2 W/kg

However, it is recommended to inform the treating physician about any implants prior to an MRI examination. Show your implant card!

(1) Zegelman M, Guenther G, Florek HJ, Orend KH, Zuehlke H, Liewald F et al. Results from the first in man german pilot study of the silver graft, a vascular graft impregnated with metallic silver. Vascular. 2009;17(4):190–6. DOI: 10,2310/6670.2009.00034.

(2) Zegelman M, Guenther G, Waliszewski M, Pukacki F, Stanisic MG, Piquet P et al. Results from the International Silver Graft Registry for high-risk patients treated with a metallic-silver impregnated vascular graft. Vascular. 2013;21(3):137–47. DOI: 10,1177/1708538113478773.

(3) Ueberrueck T, Zippel R, Tautenhahn J, Gastinger I, Lippert H, Wahlers T. Vascular graft infections: in vitro and in vivo investigations of a new vascular graft with longterm protection. J Biomed Mater Res B Appl Biomater. 2005;74(1):601-7.

(4) Ueberrueck T, Meyer L, Zippel R, Nestler G, Wahlers T, Gastinger I. Healing Characteristics of a New Silver-Coated, Gelatine Impregnated Vascular Prosthesis in the Porcine Model. Zentralbl Chir. 2005;130:71–6.

(5) Li, Y., Wu, Z. Prosthetic Vascular Grafts: Past, Present, and Future. In: Tang, YD., Yang, J., Zheng, Y., Li, Y., Zeng, Y. (eds) Application of Biomaterials in the Treatment of Cardiovascular Diseases. 2024. Springer, Singapore.

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(1) Zegelman M, Guenther G, Florek HJ, Orend KH, Zuehlke H, Liewald F et al. Results from the first in man german pilot study of the silver graft, a vascular graft impregnated with metallic silver. Vascular. 2009;17(4):190–6. DOI: 10,2310/6670.2009.00034.
(2) Zegelman M, Guenther G, Waliszewski M, Pukacki F, Stanisic MG, Piquet P et al. Results from the International Silver Graft Registry for high-risk patients treated with a metallic-silver impregnated vascular graft. Vascular. 2013;21(3):137–47. DOI: 10,1177/1708538113478773.
(3) Ueberrueck T, Zippel R, Tautenhahn J, Gastinger I, Lippert H, Wahlers T. Vascular graft infections: in vitro and in vivo investigations of a new vascular graft with longterm protection. J Biomed Mater Res B Appl Biomater. 2005;74(1):601-7.
(4) Ueberrueck T, Meyer L, Zippel R, Nestler G, Wahlers T, Gastinger I. Healing Characteristics of a New Silver-Coated, Gelatine Impregnated Vascular Prosthesis in the Porcine Model. Zentralbl Chir. 2005;130:71–6.
(5) Li, Y., Wu, Z. Prosthetic Vascular Grafts: Past, Present, and Future. In: Tang, YD., Yang, J., Zheng, Y., Li, Y., Zeng, Y. (eds) Application of Biomaterials in the Treatment of Cardiovascular Diseases. 2024. Springer, Singapore.

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The main product trademark "Aesculap" and the product trademarks "Prospace", "S4" and "activL®" are registered trademarks of Aesculap AG.
Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany